

OMRON BF511 BODY COMPOSITION MONITOR



Product Code	27293
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No
RDM (NSIS)	2556668
CND	Z12099001
EAN/UPC	4015672104051
PARAF	903604268

Description

OMRON BF511 BODY COMPOSITION MONITOR

OMRON BF511 provides full body measurement with clinical accuracy, it is validated and classified as a medical device, even for kids.

It offers a comprehensive set of readings including body and visceral fat, skeletal muscle level, BMI and resting metabolism.

Measurement with 8 sensors on both hands and feet for complete body measurement.

It uses the Bioelectric Impedance (BI) method, fat is shown as a percentage of body weight.

Multilanguage package and manual (GB, IT, FR, DE, NL, RU, Arabic).

Battery operated: 4xAA (included).