

SECA 745 BABY SCALE Class III - mechanical - 16 kg



Product Code	27294
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	76846
CND	Z12080402
EAN/UPC	4012030013550
PARAF	970930222
GMDN	173958

Description

The sturdy classic: with its softly curved weighing tray, the stable all-steel casing and the precise, smooth running weights, this model is the indestructible original baby scale.

- Metrological Approval for use in medical field
- Calibration according to class III sensitivity required by 90/384/EEC

User manual: GB, FR, IT, ES, OE, PT, PL.

Technical Specifications

- Capacity: 16 kg
- Sensitivity: 10 g
- Size: 550 x 180 x h 290 mm
- Weight: 6.4 kg
- Made in Asia