

SECA 725 BABY SCALE - mechanical - 16 kg



Product Code	27299
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No

RDM (NSIS)	1403040
CND	Y0399
EAN/UPC	4012030013567

Description

The sturdy classic: with its softly curved weighing tray, the stable all-steel casing and the precise, smooth running weights, this model is the indestructible original baby scale.

User manual: GB, FR, IT, ES, OE, PT, PL.

Technical Specifications

- Capacity: 16 kg
- Sensitivity: 5 g
- Size: 550 x 180 x h 290 mm
- Weight: 6.4 kg
- Made in Asia