

SIRIO SCALE with height meter - 150 kg



Product Code	27305
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC	8023279273052
---------	---------------

Description

SIRIO SCALE - 150 kg mechanical with height meter (85-210 cm)

Aluminium rod, height range 85 to 210 cm (division 1 mm)

Sturdy and reliable dial scale to warrant exact and precise measurements.

Large dial for easy reading.

User manual in GB, FR, IT, ES, PT, DE, GR, Arabic.

Technical Specifications

Capacity: 150 kg

Sensitivity: 500 g

Weight: 13 kg

Dimensions: 68 x 27.5 x 108 cm

Platform: 37 x 27 cm

Package size: 106 x 34 x 29 cm