

SOEHNLE 8352 MULTINA DIGITAL BABY AND CHILD SCALE



Product Code	27306
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279273069
PARAF	903170468

Description

SOEHNLE 8352 MULTINA DIGITAL BABY AND CHILD SCALE

Well-designed, easy to handle plastic baby scale.

- suitable for children up to 20 kg removing baby tray
- easy to read LCD-display with 13 mm figures
- simple operation with one button
- tare function to calibrate for baby blankets/textiles
- fully-electronic scale: no wear and tear parts, no maintenance
- hold function: weight value will be held for 10 seconds after unloading
- automatic power-off

Made in Germany.

Manual in: DE, GB, FR, IT, ES, NL, DK, SE, HU, FI, PT.

Technical Specifications

Capacity: 20 kg
Sensitivity: 10 g
Dimensions: 601x346x92 mm
Weight: 3.5 kg
Battery: 9 V
Warranty: 2 years