



Product Code	27313
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279273137
PARAF	903118053

Description

GIMA ELECTRONIC BABY SCALE

Handy, reliable and well-designed electronic baby scale.

Provided with large LCD display 66x28 mm and kg/lb switch.

Supplied with height tape meter and washable sheet.

Multilingual user manual and gift box: GB, FR, IT, ES, PT, DE, Arabic

LCD digit height: 18 mm

Technical Specifications

Capacity: 20 kg / 44 lbs
Sensitivity: 10 g
Dimensions: 595x340x90 mm
Weight: 4 kg
Battery operated: 4x1.5 V AA
Auto-off function