

SECA 416 - INFANTOMETER



Product Code	27316
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No
RDM (NSIS)	93686
CND	V0399
EAN/UPC	4012030007092
PARAF	903118103

Description

SECA 416 - INFANTOMETER

Infantometer for measuring babies up to two years old. Sturdy construction and stability makes it quick and easy to use.

Rounded sides guide the baby easily and safely into the right position.

Permanently mounted with rollers on two guide rails, the footpiece slides smoothly along the scale.

Footpiece is lockable so the baby can be removed and the result noted later.

Measuring range: 33-100 cm. Graduation: 1 mm.

The unique V-shape automatically places the baby in the right measuring position.

Size: 110x16.9x40.2 cm.

Weight: 3.8 kg.