

SENSOR PAD for Body Fat Analyzer



Product Code	27319
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
RDM (NSIS)	2462856
CND	N010299
EAN/UPC	8033695722010
GMDN	61020

Description

Gelled electrodes FIAB have been carried out for electrotherapy applications.
Indicated in fat analysis.
The highly flexible material allows an easy positioning of electrodes.
The solid electrode gel preserves its adhesivity and allows repositioning on the same patient.
Electrodes can be easily repositioned after use, thanks to their PVC silicone liner.

Latex-free products.