


PLASTIC SECURITY SEAL WITH STEAM INDICATOR - blue



Product Code	37332
--------------	-------

Unit of sale	box of 100 pcs
--------------	----------------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279373325
---------	---------------

Description

PLASTIC SECURITY SEALS - 2.7 cm LONG PIN

Security seal with steam indicator - blue

Container seals for steam sterilization processes.
Conform to ISO 11140-Type 1.