

SOEHNLE ELECTRONIC HEIGHT ROD



Product Code	27333
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	4006484500551

Description

SOEHNLE ELECTRONIC HEIGHT ROD

Measuring range: 50-240 cm (20-94.5 inches), division: 1 cm.

Manual in: DE, GB, FR, ES, IT, GR, CZ.

LCD digit height: 9 mm
Size: 345x92x28 mm
Power supply: 2x1.5 V AA batteries
Weight: 300 g