

SECA 203 CIRCUMFERENCE MEASURING TAPE



Product Code	27345
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I M

UK-REP	No
CH-REP	No
RDM (NSIS)	93583
CND	V0399
EAN/UPC	4012030007238

Description

SECA 203 CIRCUMFERENCE MEASURING TAPE with extra Waist-To-Hip-Ratio (WHR) calculator

Retractable and ergonomic measuring tape for accurate determination of body circumference. Measuring range: 15-205 cm, division 1 mm. High-quality metal scale to obtain the Waist-To-Hip-Ratio (WHR).
The WHR is an indirect method of determining the amount of abdominal fat tissue.
The ratio allows conclusions about the distribution of fat stores and thus contributes to a diagnosis.

Medical device, Class I.
Multilingual manual: GB, FR, IT, ES, DE, PT, PL.

Diameter: 70 mm
Height: 28 mm
Weight: 75 g

Display to determine body fat distribution.