

# PLUMBING WEIGHT with wooden case



Product Code	27355
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

EAN/UPC	8023279273557
---------	---------------

GTIN	08023279273557
------	----------------

## Description

### PLUMBING WEIGHT

Made of stainless steel (magnet) with wire length 150 cm.

Supplied in wooden case.

Size: 6.2x2.3 cm.

Weight: 140 g.