



Product Code	27373
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	102604
CND	Z110702
EAN/UPC	8023279273731
PARAF	903112858
GMDN	36488

## Description

GIMA light boxes are designed according to clinical standards for an extra white light for an accurate and precise reading of any detail.

Made in Italy.

- External dimensions
  - Length: 43 cm
  - Height: 127 cm
  - Depth: 12 cm
- Weight: 14 kg
- 2 vertical panels

Incorporates a 3 mm acrylic screen compliant with safety norms (fuses, Schuko plug, green light On/Off, ...).

Can be used on the table placed horizontally, vertically, depending on model or wall mounted.

Body made of oven-fired epoxy powder coated aluminium and injected ABS (colour RAL 7035).

## Technical Specifications

Colour temperature: Cool white 6,500 °K  
 Power supply: 230 V - 50 Hz (60 Hz available with surcharge)  
 Protection fuse: 2x1A  
 Bipolar light switch  
 Schuko plug with 2 m cable  
 All components UL, DVE, IMQ approved