

"ARMAGNAC" CHART - 5 m - not illuminated



Product Code	27390
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	102748
CND	Z12120185
EAN/UPC	8023279273908
PARAF	903120537
GMDN	16800

Description

Armagnac Optometric chart

- Not illuminated model
- Reading distance: 5 meter

GIMA optotypes have been designed following clinical recommendations to give an accurate and precise reading of any detail.
Manufactured with acrylic screen.

Made in Italy.

Technical Specifications

- Not illuminated model
- Reading distance: 5 meter
- Viewing area: 24 x 61 cm
- External dimensions: 29 x 66 x 2.5 cm
- Weight: 2 Kg