

PEDIATRIC BED 1-4 years



Product Code	27399
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2537673
CND	V080602
EAN/UPC	8023279273991
GTIN	08023279273991
GMDN	37010

Description

PAEDIATRIC BED 1-4 years - ivory

Epoxy powder painted mainframe made of steel tube Ø 30 mm, ready to be hooked into the mattress platform.
Mesh mattress platform made of oval tube 40x20 mm, fixed to the mainframe with removable bolts.
Side rails 62 cm h. sliding on chromed steel guides, with safety device.
Legs with 4 castors Ø 80 mm, two with brake.

Weight 35 kg.
External measurement 140x70xh 130 cm
Inner size 130x65xh 63 cm.

Made in Italy.