

EXAMINATION COUCH - white painted, grey mattress



Product Code	27402
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2537687
CND	V080603
EAN/UPC	8023279274028
PARAF	903605727
GTIN	08023279274028
GMDN	38458

Description

EXAMINATION COUCH - white painted - grey

Suitable for medical office or small surgeries.
Easy to assemble, following enclosed instructions (GB, FR, IT, ES, PT, DE, SE).
Comes complete with mattress.

Upholstery in two sections: one 131x62xh 6 cm fixed and a 51x62xh 6 cm backrest adjustable through different positions.
Load: 120 kg.
Net weight 26.4 kg
Delivered in kit form (173x16x74 cm).

Backrest adjustable 0° to 55° through different positions

Fireproof, simulated leather, washable upholstery in two sections

Made in Italy.