

ALUMINIUM SUITCASE COUCH



Product Code	27406
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2537691
CND	V080603
EAN/UPC	8023279274066
GTIN	08023279274066
GMDN	38458

Description

ALUMINIUM SUITCASE COUCH with hole

Light, foldable, two sections, suitcase couch with aluminium frame.
Central double legs for improved safety and block system on the hinges.
Head rest with mouth/nose hole (hole 23x12 cm).
Upholstered and covered with simulated leather.
Max load: 100 kg.
Weight: 13 kg.
Packing: 94x21.5x66.5 cm

Manual: GB, FR, IT, ES, PT, DE, GR

Size: 183x63xh 76 cm

Made in Italy.