

# SUITCASE COUCH



Product Code	27407
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2537691
CND	V080603
EAN/UPC	8023279274073
PARAF	903112922
GTIN	08023279274073
GMDN	38458

## Description

### SUITCASE COUCH - black

Chromed steel and folding legs.  
Upholstered covering in washable material.  
Resin padded.  
Max load: 100 kg.  
Weight: 16.5 kg  
Supplied with a velvet bag for easy transport.  
Packing: 95x18x73.5 cm  
Manual: GB, IT, NL.

Couch size:  
187.5x64xh 71 cm (open)  
93x68xh 14 cm (folded)

Mattress thickness 3 cm

Made in Italy.