

# WOODEN COUCH WITH HOLE



Product Code	27417
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279274172

## Description

WOODEN COUCH - with face hole - beech

MAXIMUM COUCH LOAD: 200 KG

Comfortable couches, ideal for physical therapy, esthetics and wellness centres.

Superior quality wood.

Back rest adjustable in 13 positions.

Resilient foam padding covered with washable skay.

Padding in beige colour.

Size: 190x70xh 75 cm.

Delivered in kit form. Easy assembly, tool provided.

Mouth nose hole plug included

Made in Italy.