

# COUCH ROLL HOLDER for round or squared 30mm metal leg



Product Code	27418
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2275292
CND	V9099
EAN/UPC	8023279274189
PARAF	903112959
GTIN	08023279274189
GMDN	38458

## Description

Adjustable roll holder suitable for round ( $\varnothing$  27-40 mm) or square (25-25 mm x 35-35 mm) metal couch leg.