

# MICROEMBOSSED GLUED 2 PLIES PEDIATRIC COUCH ROLL - 46m x 50cm



Product Code	27423
Unit of sale	9 rolls
Minimum order	1
Type	No medical device
Class	

EAN/UPC 8017414001327

## Description

MICROEMBOSSED GLUED 2 PLIES PEDIATRIC COUCH ROLL - 46m x 50cm - box of 9

Decorated medical rolls for pediatric use. Produced in pure cellulose, dermatologically tested.  
Hypoallergenic inks.  
Each box contains 9 rolls with 5 assorted themes.

Material: Microembossed glued  
Colour: white+decoration  
Ply: 2 ply  
Weight: 2x17.5 g/m<sup>2</sup>  
Length: 46 m  
Height: 50 cm  
Number of sheets: 135  
Perforation: 34 cm  
Rolls per box: 9