

# GIMA 3 TROLLEY without guard-rail - medium



Product Code	27426
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279274264
GTIN	08023279274264

## Description

GIMA 3 TROLLEY - without guard rail - medium 70x50xh 78 - weight 15 kg.

Practical trolleys with aluminium body and handles.  
Two white bilaminated shelves 18 mm thick, black PVC bordered, without guard rail.  
Two large lateral handles facilitate transport or lift up of the trolley.  
4 revolving nylon castors Ø 80 mm.  
Max load (x shelf): 40 kg.

Delivered in kit form.

Made in Italy.