

GIMA 2 TROLLEY with guard-rail - small



Product Code	27429
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279274295
PARAF	903605816
GTIN	08023279274295

Description

GIMA 2 TROLLEY - small: 60x40xh 77 cm - weight 11 kg

Practical trolleys with aluminium body and handles.
Two white bilaminated shelves 18 mm thick, black PVC bordered, with guard-rail.
Two large lateral handles facilitate transport or lift up of the trolley.
4 revolving nylon castors Ø 80 mm.
Max load (x shelf): 40 kg.

Delivered in kit form.

Made in Italy.