



Product Code	27430
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279274301
PARAF	903112961
GTIN	08023279274301

Description

INOX TROLLEY

AISI 304 Stainless steel trolley with plastic drawer internally divided in 4 sectors. 2 shelves with upstands (3,5 cm) on four sides. Revolving rubber castors Ø 80 mm (two with brake).

Size: 70x50xh 82 cm.
 Max load (x shelf): 25 kg.
 Delivered in kit form.

Made in Italy.