



Product Code	27436
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279274363
PARAF	903605879
GTIN	08023279274363

Description

DRESSING TROLLEY

Dressing trolley made of AISI 304 stainless steel with 2 shelves (70x50 cm), 1 drawer (35x46xh 12 cm), basin holder and basin (Ø 32xh 9 cm).

4 revolving castors Ø 100 mm.

Size: 70x50xh 85 cm - weight: 18 kg.

Max load (x shelf): 25 kg.

Delivered in kit form.

Made in Italy.