



Product Code	27439
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8023279274394

Description

REX TROLLEY

Complete stainless steel dressing trolley.
Supplied with 3 shelves, waste bin, drawer, support for bottles, guard rails on the two bottom shelves, 4 wheels, Ø 100 mm (2 with brakes).

Load 50 kg.
Size: 79x50xh 120 cm (handle included).
Drawer size 64x43xh 9.5 cm.
Delivered in kit form (114x54x37 cm).

Made in Italy.