



Product Code	27442
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No

EAN/UPC	8023279274424
GTIN	08023279274424

Description

INSTRUMENT TROLLEY

Made of satined welded AISI 304 stainless steel. With 4 sliding drawers (internal, 49x33xh 18 cm) Top drawer lockable, upper shelf with safety edges to carry instruments. 4 revolving castors Ø 80 mm

Size: 60x40xh 110 cm

Weight: 37 kg

Made in Italy.