

FIRENZE WOODEN COUCH - wenge



Product Code	27450
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279274509

Description

FIRENZE WOODEN COUCH - wenge

MAXIMUM COUCH LOAD: 250 KG

Comfortable couches, ideal for physical therapy, esthetics and wellness centres.
Superior quality wood.
Back rest adjustable in 13 positions.
Resilient foam padding covered with washable skay.
Padding in beige colour.

Size: 190x70xh 75 cm.

Delivered in kit form. Easy assembly, tool provided.

Couch roll holder included can be hanged on back or front.

Made in Italy.