

DELUXE TROLLEY with drawer 58 x 40 X H 19.5 cm



| | |
|---------------|-------------------|
| Product Code | 27461 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | No medical device |
| Class | |

| | |
|--------|----|
| UK-REP | No |
| CH-REP | No |

EAN/UPC 8023279274615

Description

DELUXE TROLLEY with drawer

Size: 65x45xh 80 cm

Drawer size: ext. 58x40xh 10.5 cm, int. 51.5x38.5xh 8.5 cm.

Top quality trolleys made of lightweight anodized special alloy with detachable stainless steel shelves.

4 revolving rubber castors Ø 80 mm.

Max load (x shelf): 20 kg.

Instructions: GB, FR, IT, ES PT, DE, GR.

Delivered in kit form.

Made in Italy.