

LAUNDRY TROLLEY - stainless steel



Product Code	27465
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279274653
PARAF	987476621

Description

LAUNDRY TROLLEY - stainless steel

Trolleys to collect dirty laundry while distributing clean linen to the patient. Stainless steel or painted steel frame. Equipped with 4 Ø 125 mm castors (2 with brake) and 4 rubber bumpers at corners.

Size: 135x65xh 102 cm.

Wrapping bag not included

Made in Italy.