

# OVERBED TABLE - moulded tray



|               |                   |
|---------------|-------------------|
| Product Code  | 27469             |
| Unit of sale  | 1 pc              |
| Minimum order | 1                 |
| Type          | No medical device |
| Class         |                   |

|         |               |
|---------|---------------|
| UK-REP  | Not required  |
| CH-REP  | Not required  |
| EAN/UPC | 8023279274691 |

## Description

OVERBED TABLE - moulded tray

Frame made of chrome plated steel tube.  
Serving tray made of thermoformed anti-shock plastic material with retaining edges and pre-formed bottle, glass and cutlery holder.

Tray dimensions: 58 x 39 x h 5 cm  
Adjustable height: 75 to 114 cm  
Base with 4 Ø 50 mm castors.  
Delivered in kit form. Carton size: 59 x 22 x 63 cm.

Made in Italy.