

OVERBED DIFFUSION TABLE - white



Product Code	27475
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2044511
CND	Y1803
EAN/UPC	3701127706130

Description

OVERBED DIFFUSION TABLE - white

Overbed table with high quality, laminated, fire resistant tabletop with an edging strip.
Epoxy coated "C" steel frame allows for closer access in a chair.
Defivered ready to use and easy to fold flat (15 cm) for ease of storage.

Weight: 8 Kg
Max load: 10kg

Table top: 60 x 40 cm
Adjustable height: from 61 to 94 cm
Four castors, two with brake.

Made in France.