

DRESSING TROLLEY - medium



Product Code	27484
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279274844
PARAF	903607012

Description

DRESSING TROLLEY - Medium
Colour: yellow
Size of trolley (without accessories): 102 x 64 x h 100 cm
Number of drawers: 8
Size of drawers:
- 45 x 48 x h 13 cm (3)
- 45 x 48 x h 27 cm (1)
- 30 x 48 x h 13 cm (3)
- 30 x 48 x h 27 cm (1)

These multifunctional professional trolleys have painted steel structure and external plastic finishing.
Provided with 4, 8 or 9 drawers, one lateral waste container (manual opening, capacity 15 l) and one extra lateral table.
Four swivelling rubber castors Ø 125 mm (2 with brake).
Equipped with instrument tray and upper dispenser (10 tilting bins).

Made in Italy.