

ELITE OVERBED TABLE



Product Code	27487
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

EAN/UPC	8023279274875
---------	---------------

Description

ELITE OVERBED TABLE

High quality overbed table with crank adjustment system.
It features a base that allows closer access.
It is equipped with 4 hidden castors to easily maneuver it anywhere.
It can enhance patient well-being and comfort to elevate the quality of healthcare.

Size: 800x400xh 615-905 mm.
Max. load: 25 kg.
Weight: 11.5 kg.
Packaging size: 820x90x430 mm

Delivered in kit form.

MDF (medium density fiberboard) table top 800x400 mm
Height adjustment (615-905 mm)
Steel frame
Folding handle
Castors Ø 16 mm

Multilingual manual: GB, FR, IT, ES, PT, DE.