

3-SECTION MULTIFUNCTIONAL BED - blue



Product Code	27502
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2254540
CND	V0899
EAN/UPC	8055136920045

Description

3-SECTION MULTIFUNCTIONAL BED - blue

Strong and stylish structure in oval tube 50x25 mm, painted with epoxy powder, metallic grey colour. The three sections top has high density padding, covered with washable, fire proof, leatherette material Class 1 IM. Back-rest and leg-rest section are fully adjustable by self mechanical sytem.

Maximum load: 150 kg - Weight: 36 kg.

Seat: 62x60 cm

Size: 185 x 60 x h 85 cm.

Delivered in kit form.

Made in Italy.