



|               |                |
|---------------|----------------|
| Product Code  | 27515          |
| Unit of sale  | 1 pc           |
| Minimum order | 1              |
| Type          | Medical device |
| Class         | I              |

|            |               |
|------------|---------------|
| RDM (NSIS) | 2248544       |
| CND        | Z12011201     |
| EAN/UPC    | 8023279275155 |

Description

Gynex professional chair offers fuctionality, aesthetics, ergonomics.

It can satisfy all patient and user needs for to a professional gynaecologic visit. Gynex professional is provided with versatile accessories: back paper-roll dispenser, metal basin, chromed handles, leg-supports.

3 engines allows, in addition to lifting, electric backrest and legrest movements.

The combined movement of backrest and seat by means of self-balancing system (operated by a silent low tension actuator of 2500 N) grants maximum comfort to the patient.

The padded and pelted armrests allow a relaxing protecting sitting position, armrest are individually revolving and removable. Backrest is provided with removable height-adjustable headrest. The wide and soft backrest has independent adjustment movement through silent low tension actuators (1000 N each one).

This chair is moved by a silent low tension engine which allows a vertical excursion of 40 cm from zero position.

The electrical drives are controlled by an ergonomic double switch control, connected to the chair by a cable. Gynex professional can be easily transformed into an examination bed (horizontal position).

Made in Italy.

Technical Specifications

- Power supply tension: 230 - 50 hz
- Motor tension: 24 V
- Lifting motor push: 6000 N
- Backrest motor push: 2500 N
- Legrest motor push: 2500 N
- Lifting capacity: 180 kg
- Maxium power: 90 VA
- Fuser protection: 1.6 A

Standard accessories

- Padded leg-rest (colour coordinated)
- Paper-roll holder
- Metal basin
- Remote control