

GYNEX BED CHAIR with roll holder - blue



Product Code	27523
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2248539
CND	Z12011201
EAN/UPC	8023279275230

Description

GYNEX - GYNAECOLOGICAL CHAIR

GYNEX: a chair offering a unique combination of functional, aesthetic and working characteristic. Thanks to its self-balancing seat and backrest, GYNEX can be adjusted to any build for the greatest patient comfort. The chair is provided with an electromechanical lifting device, having a power supply at 24 V. The electric actuators enable the movement up and down from 57 cm off the floor ("0" position) up to a maximum of 96 cm. The backrest is supplied with adjustable headrest and paper roll holder. The legrest can be independently adjusted and it has an automatic-locking device. The chair can be flattened to a bed position for examination activity, the armrest can be removed and independently oriented. The electric control of the chair for upright position is operated by means of a foot pedal.

Made in Italy.

Technical Specifications

- Power supply tension: 230V - 50 hz
- Motor tension: 24 V
- Lifting motor push: 6000 N
- Backrest push: 800 N (gas spring)
- Lifting capacity: 180 kg
- Maxium power: 90 VA
- Fuser protection: 1.6 A

Standard accessories

- Poliuretane leg-rest
- Paper-roll holder
- Basin (plastic)