

LED LIGHT SOURCE



Product Code	27544
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2208042
CND	Z12020402
EAN/UPC	8023279275445

Description

LED LIGHT SOURCE

Cold white light (5,500 °K, 200 lumen) for endoscopy use. It can be connected to endoscopes without need of F.O. cable. Handle connect directly to rigid endoscopes (Storz connection). 3 Watt power LED.

Size: 30x25x12 cm
LCD display: 12.8x6.4 cm
Adjustable light

Made in Italy.