



Product Code	27561
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	1229071
CND	Z129008
EAN/UPC	8023279275612

Description

DONOR WHEELCHAIR

Donor wheelchair with structure made of epoxy painted oval and round steel tube.
 Backrest adjustable by gas pump.
 Padded leg rest manually adjustable.
 Metal footrest extractable. Armrest with articulations, manually adjustable in height.
 Padding made of high-density polyurethane foam, fireproof Class 1 IM.
 Covering made of washable, plastic material (leatherette), fireproof Class 1 IM.
 Back push handle.
 Mattress overall size: 166x55 cm.
 Carton size: 185x67xh 58 cm
 Four wheels Ø 100 mm, 2 directional and 2 with brake.
 Net weight: 53 kg.
 Gross weight: 55 kg.
 Maximum load: 150 kg.

Supplied with blanket tray with belts to fix oxygen bottle, IV stand and chart holder support.

Made in Italy.