



Product Code	27567
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2255152
CND	V080202
EAN/UPC	8055136920113

Description

DONOR ARMCHAIR WITH WHEELS - white

Strong and stylish armchair with structure made of oval, round and square metal tubing fully coated with white epoxy powder (legs oval tubing 50x25x1.5 cm; central frame square tubing 30x30x1.5 cm, cushions round tubing 20xØ 1.5 cm). High density polyurethane foam platform with wooden base covered by leatherette upholstery. Flame resistant materials Class 1 IM.

Synchronized adjustment of the two cushions by a gas spring (cylinder/piston) allows a quick moving into TR/RTR position.

Backrest inclination from +60° to -5° in Trendelenburg position, leg rest from -80° to +5°.

Reclinable arm-rests can be height adjusted or rotated up to 180° when removed.

4 castors, 2 with directional brake ensure mobility.

Load: 150 kg.

Size: 175x55xh 51 cm.

Multilingual manual: GB, FR, IT, ES, PT, DE, SE, PL, RO, GR

Made in Italy.