

HYPODERMIC NEEDLE 23G 0.6x25 mm - sterile



Product Code	23758
Unit of sale	box of 100 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2357576
CND	A0101010102
EAN/UPC	8032764151157

Description

STANDARD HYPODERMIC NEEDLES 23G - STERILE

Standard hypodermic needles, made in stainless steel AISI 304 with triple bevelling and totally non-toxic and non-pyrogenic. Each element of the needle is treated for painless injection and ensuring minimal trauma during penetration.

- Size: 23G - 0.6x25 mm
- Reference colour: turquoise

- Box of 100 pieces
- Individually packaged in blister
- Box and instructions: GB, IT, ES, PT.