

THER TRENDELENBURG TABLE with armrest - electric - blue



Product Code	27599
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

RDM (NSIS)	2346891
CND	V080603
EAN/UPC	8023279275995

Description

THER TRENDELENBURG TABLE WITH ARMREST - electric - blue colour

Examination and treatment tables with electric height adjustment by a foot switch (low tension 24 V).

Features: three mobile sections adjustable by gas spring. The table can easily assume the “flexion”, “relax” and “trendelenburg” position even if the patient is lying on the table. The large head-section (with mouth-nose hole) can be adjusted positively and negatively.

Chair position.

Safe working load 200 kg.

Power 230 V - 50/60 Hz.

Supplied with mouth-nose hole plug.

Made in Italy.

Technical Specifications

- Lifting capacity: 200 kg
- Total Length: 200 cm
- Table width: 68 cm.
- Height: adjustable from 47 to 97 cm
- Adjustable head-section: +85°-35°