

# BED 1 ARTICULATION - with wheels 100 mm



Product Code	27601
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2537730
CND	V080602
EAN/UPC	8023279276015
PARAF	903113191
GTIN	08023279276015
GMDN	34873

## Description

### BED WITH ADJUSTABLE HEAD SECTION

- BED WITH 1 ARTICULATION with castors

Painted steel frame with net.

1 Pipe frame with net for head movement, 200 x 85 cm.

2 Heads in chromed pipe complete with MDF bilaminated coated panel.

Movement by a crank.

Protection of end run.

Rubber wheels of Ø 100 mm, with brakes.

Size: 206 x 89 x h 90 cm.

Height of frame from bottom 52/55 cm.

Size of frame: 200 x 85 cm.

Weight 32 kg.

Load: 110 Kg.

Made in Italy.