

# BED 3 ARTICULATIONS with wheels



Product Code	27611
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2537730
CND	V080602
EAN/UPC	8023279276114
PARAF	903113266
GTIN	08023279276114
GMDN	34873

## Description

### BED WITH ADJUSTABLE HEAD AND FOOT SECTION

- BED WITH 3 ARTICULATIONS with castors

Painted steel frame with net.  
3 bed frames with nets in painted steel pipe, 200 x 85 cm.  
2 Heads in chromed pipe complete with MDF bilaminated coated panel.  
Movement by 2 cranks.  
Protection of end run.  
Rubber wheels of Ø 100 mm, two of which with brake.  
Size 206 x 89 x h 90 cm.  
Height of frame from bottom 52/55 cm.  
Size of frame: 200 x 85 cm.  
Load: 110 Kg.  
Weight 42 kg.

Made in Italy.