

KING EXAMINATION COUCH - water green



Product Code	27614
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2254328
CND	V0899
EAN/UPC	8023279276145

Description

Strong and stylish structure in oval tube 50 x 25 mm painted with metallic grey epoxy powder. Top in washable, fire proof, leatherette material class 1 IM with high density padding. Head section adjustable by mechanical system, maximum 75°.

Mattress thickness: 4.5 cm
Maximum load: 150 kg.
Size: 190 x 60 x h 80 cm
Delivered in kit form.

Made in Italy.