

ELECTRIC HEIGHT ADJUSTABLE TREATMENT TABLE - beige



Product Code	27634
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	Yes

RDM (NSIS)	2537731
CND	V08060199
EAN/UPC	8023279276343
PARAF	903113304
GTIN	08023279276343
GMDN	38447

Description

GIMA ADJUSTABLE HEIGHT ELECTRIC TABLE - footswitch

- Electric table - footswitch
- Colour: beige
- Mouth-nose hole plug - optional: 27620

Sturdy and highly stable examination and treatment table with white metal painted frame.

Adjustable height by electric footswitch.

Fixed leg-section and manually adjustable head-section with mouth-nose hole.

Available in 12 colours.

Fireproof padding.

Class I Medical Devices, complies with MDR 2017/745.

Lifting capacity: 135 kg.

Length: 193 cm.

Weight: 54 kg.

Power supply 100-240 V, 50-60 Hz.

Multilingual manual: GB, FR, IT, ES, PT, DE, GR.

Made in Italy.