



Product Code	27645
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2302554
CND	V0899
PARAF	987472646
GTIN	18054117420458

## Description

### LIFTING STAND

Tubular painted steel frame with adjustable chromed steel joint.  
 Can be fully disassembled to reduce size.  
 Three main height adjustments.  
 Triangle hangs from belt with adjustable length.

Adjustable height 182 to 200 cm  
 Base 87x70 cm

## Technical Specifications

Overall length: 864 mm  
 Overall width: 700 mm  
 Overall height: 2000 mm  
 Ring center: -  
 Load\*: 80 kg  
 \*Lifting stand load considers half body lift