


S/S I.V. STAND FOR BED - 4 hooks



| | |
|---------------|----------------|
| Product Code | 27646 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | I |

| | |
|--------|-----|
| UK-REP | Yes |
| CH-REP | Yes |

| | |
|------------|----------------|
| RDM (NSIS) | 2753826 |
| CND | V9099 |
| EAN/UPC | 8023279276466 |
| PARAF | 903113367 |
| GTIN | 08023279276466 |
| GMDN | 36069 |

Description

4 Hook IV stand for bed - stainless steel.

Height from 41 to 64 cm.