



Product Code	27664
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No
RDM (NSIS)	2517662
CND	V08068002
EAN/UPC	8023279276640

Description

Lift pole made of steel frame with coating surface.
The frame is composed of 4 parts, easy to assemble.
Supplied with hand ring.

Technical Specifications

Overall length: 725 mm
Overall width: 790 mm
Overall height: 1935 mm
Ring center: 1520 mm
Load*: 75 kg
*Lifting stand load considers half body lift