

UNCUFFED TRACHEOSTOMY TUBE diam. 4 mm



Product Code	57671
Unit of sale	box of 10 pcs
Minimum order	1
Type	Medical device
Class	II A

UK-REP	No
CH-REP	No

RDM (NSIS)	2724996
CND	R0105010102
EAN/UPC	6972091921356

Description

Uncuffed tracheostomy tubes Ø 4.0 mm

Curve tubes designed to rapidly provide an airway within seconds during an emergency by inserting them into a tracheostomy stoma, a hole created in the neck and windpipe.

The hollow tube whether with or without cuff, is directly placed into the trachea through a surgical incision, using a wire-guided progressive dilation technique.

- tube is made of medical grade, non-toxic, clear and kinkresistant PVC to protect delicate mucosal tissue
- smooth rounded obturator's tip reduces trauma during insertion
- neck plate with printed size for easy identification
- standard connector
- radio-opaque line helps to determine exact tube location

Tube includes: inflation line, valve, obturator, neck plate, connector and tube ties.

Sterile, single patient. Class IIa. DEHP-free.

Made with no natural rubber or latex.